





In Omero.insight you create projects/dataset (1)
and upload data (2)



In Omero.insight Administration you create groups (1)
And define permissions (2). In Edit group (3) you add members



2




In Omero webclient you open files by selecting the image and press (1) and you have 3 options:

Image viewer

OMERO i.viewer : <https://www.openmicroscopy.org/omero/iviewer/>

OMERO.figure: <https://figure.openmicroscopy.org/>



Note! You can also open "Figure" from the top bar and select images from there

Follow this tutorials above to make annotations and figures

